COLLECTION KIT INSTRUCTIONS

INTENDED USE

This kit is intended for the collection and preservation of dried blood spot samples for the purpose of laboratory testing by Simple HealthKit laboratory.

For accurate results, please read all instructions carefully before using the kit.

IMPORTANT

Create Your Online Account

Before sample collection, you must first create an account in the patient portal.

You will then register each Sample Card to your account using both the Sample Card ID found on the Sample Card and your Provider's ID. If you do not have your Provider's ID, contact them directly.

The lab will not be able to process your sample if you skip this step.

LIMITATIONS

Do not use if you have a bleeding disorder or are taking medications that increase the risk of bleeding.

Do not use if you are at risk of fainting when providing blood samples.

Do not use provided lancets on broken, damaged, ulcerated, scarred, or infected skin.

There is a possibility that using these lancets may cause bruising, tenderness, soreness, pain, nerve damage, scarring on the skin, damage to the bone, or infection. Seek medical attention if you have any concerns.

Testing using Clara Spot™ is available with an order from your health care provider. This testing has been validated for use in individuals 18 or older.

Keep out of reach of children. This kit contains small parts that may pose a choking hazard.

Store your Sample Cards in a cool dry place between 59°F (15°C) and 77°F (25°C) away from sunlight and moisture. Do not freeze.

TO COLLECT EACH SAMPLE

Each kit includes:

ALCOHOL WIPES (2)

LANCETS (2)

SAMPLE CARD

BANDAGES (2)

SAMPLE BAG

RETURN ENVELOPE
TO RETURN ALL SAMPLES

DO NOT USE any components other than what is provided in the Collection Kit.

DO NOT USE this kit if any components are missing, damaged or expired—visit claraspot.com/support for support. The expiration date can be found on the kit package.

Note: You may be collecting multiple samples—do not discard the instructions.

01

REGISTER YOUR CARDS

Print your name, date of birth, and date of collection in the provided fields on your Sample Cards.

Your name and date of birth on the Sample Cards must match the information in the portal. Unlabeled Sample Cards may be rejected or delayed.

Scan the QR code or visit claraspot.com to log into the patient portal. If you do not have an account, click on “Create Account” and follow the steps to make your account.

Once logged in to your account, the portal will direct you to “Register Your Kit.” If it does not, click on “Register a Kit.” You must register each Sample Card using its unique Sample Card ID and your Provider ID. If you do not have your Provider ID, contact your provider directly for it.

For additional information and support, visit claraspot.com/support.

02

PREPARE AND COLLECT YOUR SAMPLE

Lay out collection materials: alcohol wipe, lancet, Sample Card, and bandage. Open and place your Sample Card on a flat surface.

Wash your hands with warm water and soap for at least 1 minute and dry completely.

Select and clean a finger for sample collection. We recommend using your ring finger (4th finger) or middle finger. Use the provided alcohol wipe to clean your fingertip. **Let your fingertip dry completely. DO NOT blow on it.**

Shake your hand vigorously for 30 seconds to stimulate blood flow.

TIPS FOR A SUCCESSFUL SAMPLE COLLECTION

To increase blood flow try the following:

Aim to drink 64 oz of water the day before, and **16 oz of water 30 to 60 minutes** prior to sample collection.

Avoid caffeine, nicotine or tobacco 2 hours before collecting sample(s).

Warm your hands prior to sample collection.

Keep your hands **below your heart** during sample collection.

COLLECT YOUR SAMPLE

Twist the top off of the lancet and discard. For the finger you will prick, place the white tip of the lancet against the right or left side of the clean fingertip. Firmly press the lancet against your finger to puncture the skin. You should see a blood drop when you remove the lancet.


Gently massage from below the puncture site toward the tip to encourage blood flow. Allow several seconds for a large blood drop to form. Hold your finger pointing down over the first spot on the Sample Card.


Let a large blood drop fall onto the center of the spot and spread to the dotted line. Add drops until the spot is filled and soaks to the back of the card. This may take 2-3 drops. Make sure you have fully filled the first spot before moving on to the next. Completely fill at least 3 spots.

DO NOT press or smear your fingertip onto the card. **DO NOT** rub or smear the blood sample.

Make sure each spot on the Sample Card is completely filled, allowing enough blood to soak through to the back of the card before moving on to the next spot.

ACCEPTABLE
At least 3 completely filled spots

UNACCEPTABLE
Blood spots are overlapping

UNACCEPTABLE
Not enough blood on each spot

Note: If you cannot fill at least 3 spots with one finger, use the second lancet on a different finger. Be sure to follow the same cleaning procedure before collecting blood from the new finger.

When done collecting all blood spots, bandage your fingers if needed.

DRY YOUR SAMPLE CARDS

Allow the blood spots to air-dry completely.

Leave the Sample Card open and let the blood spots dry for **at least 20 minutes** before closing. Keep the card in a cool, dry place away from moisture and direct sunlight while it dries.

DO NOT close the Sample Card until the blood is fully dry.

Once the blood spots are completely dry, close the Sample Card by tucking the cover into the flap. Place each Sample Card into its own Sample Bag and seal the bag completely.

Note: If your health care provider asked you to collect a second sample, store the first Sample Card inside its individual Sample Bag in a cool, dry place away from moisture until you are ready to send BOTH Sample Cards to the lab.

RETURN SHIPPING CHECKLIST

- ☐ I registered my Sample Card with my Provider's ID, Sample Card ID, and the date of sample collection.
- ☐ My registered name, birth date, and the date of sample collection are on each Sample Card. This information matches the information in the portal.
- ☐ I collected the sample(s) according to the instructions.
- ☐ I filled at least 3 spots on each Sample Card.
- ☐ I placed each Sample Card in its own Sample Bag.
- ☐ I will drop off the package at a FedEx® location, or arrange for the package to be picked up by FedEx® within 24 hours of collection.
- ☐ My FedEx® tracking number is _____

03

SHIP YOUR KIT

Make sure each completed Sample Card is in its own sealed Sample Bag. Place the Sample Bag(s) into the return envelope and seal.

IMPORTANT

Drop your package off at a FedEx® location or schedule a pick up as soon as possible.

Do not use a FedEx® dropbox if outside temperatures are below 35°F (2°C) or above 90°F (32°C).

Find the nearest FedEx® drop-off location and collection hours at local.fedex.com. You can get free pick up Monday through Friday at most locations. Contact FedEx® at +1 (800) 463-3339 or on the web at fedex.com/en-us/shipping/schedule-manage-pickups.html at least one day in advance.

Keep your **tracking number** to help you track the package.

04

GET YOUR RESULTS

You'll receive an email and text notification when your results are ready. You can also log in to your account to view updates. Results are typically available 2 to 4 days after your samples are received by the lab.

Your health care provider will help you interpret results.

CLARA SPOT™ SUPPORT

Questions about your sample collection or Clara Spot™ Testing Service? Contact support at claraspot.com/support.

Clara Spot™ Kit manufactured and tested by Simple HealthKit, Inc.
446 S. Hillview Dr., Milpitas, CA 95035
Use by prescription only.
GBP-CLB-501 07/2026

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Kit de Clara Spot™ fabricado y probado por Simple HealthKit, Inc.